



ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY
1	Cart Top.....	101-013	1	12	3-Way Toggle Valve.....	120-007	2	23	Screw, 10-32 x 3/8".....	001-070	1
2	Water Relay.....	123-000	1	13	Screw, 4 x 3/8.....	001-135	4	24	Handpiece Adjustment Block.....	113-021	1
3	Routing Valve.....	014-007	1	14	Universal Handpiece Holder.....	141-006	2	25	Ferrule, 1/4".....	117-008	1
4	Plug, 1/2.....	010-013	2	15	Automatic Handpiece Holder.....	122-049	1	26	Barb, 10-32 x 1/8".....	022-010	1
5	Chassis.....	101-014	1	16	Screw, 4 x 3/4".....	001-063	2	27	Plastic Washer, 1/4".....	003-001	1
6	Plastic Elbow, 1/16.....	022-039	1	17	Screw, 6 x 1/2".....	001-072	4	28	Barb, 10-32 x 1/8".....	022-009	1
7	Gauge.....	026-003	1	18	Needle Flow Valve.....	121-000	1	29	Ferrule, 1/8".....	117-007	1
8	Nut, 6-32.....	002-004	2	19	Foot Control.....	126-000	1	30	Shuttle Valve.....	014-012	1
9	Washer, Flat, #6.....	003-032	4	20	Cart Frame.....	108-016	1	31	Plastic Tee, 1/16".....	022-049	1
10	Faceplate.....	101-015	1	21	Screw, 1/4-20 x 7/8".....	001-011	4	32	Plastic Tee, 1/8".....	022-033	1
11	Screw, 6-32 x 1/4.....	001-133	4	22	Washer, 5/16".....	003-033	4	33	Syringe.....	113-400	1

**AUTOMATIC DOCTORS CART
A-2100**

Rev 2 4/4/10