



ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY
1	Chassis.....	101-014	1	11	Universal Handpiece Holder.....	141-006	2	20	Plastic Tee, 1/8".....	022-049	1
2	Plug, 1/2".....	010-013	2	12	Automatic Handpiece Holder.....	122-049	1	21	Plastic Tee, 1/16".....	022-033	1
3	Plastic Elbow, 1/16".....	022-039	1	13	Screw, 4" x 3/4".....	001-063	2	22	Syringe Block.....	113-021	1
4	Gauge.....	026-003	1	14	Ferrule, 1/4".....	117-008	1	23	Water Relay.....	123-000	1
5	Nut, 6-32".....	002-004	2	15	Barb, 10-32 x 1/16".....	022-010	1	24	Routing Valve.....	014-007	1
6	Washer, Flat, #6.....	003-032	4	16	Plastic Washer, 1/4".....	003-001	1	25	Needle Flow Valve.....	121-000	1
7	Faceplate.....	101-015	1	17	Barb, 10-32 x 1/16".....	022-009	1	26	Screw, 10-32 x 3/8".....	001-072	1
8	Screw, 6-32 x 1/4".....	001-133	4	18	Ferrule, 1/8".....	117-007	1	27	Syringe.....	113-400	1
9	3-Way Toggle Valve.....	120-002	2	19	Shuttle Valve.....	014-012	1	28	Optional Control Head Cover.....	101-016	1
10	Screw, 4" x 3/8".....	001-135	2								

**AUTOMATIC HANDPIECE CONTROL
2100**