



ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY
1	Mounting Bracket	131-034	1	15	Screw, 4 x 3/8"	001-135	4
2	Washer, #10	003-024	2	16	Faceplate Extension	101-024	2
3	Screw, #6 x 1/2"	001-072	2	17	Screw, 4 x 1/2"	001-063	4
4	Cannister Mounting Bracket	131-038	1	18	Faceplate	101-015	1
5	Screw, 8-32 x 3/8"	001-003	2	19	Vacuum Holder	141-038	1
6	Vacuum Cannister	131-035	1	20	3-Way Toggle Valve	120-008	2
7	Cart Top	101-013	1	21	Universal Handpiece Holder	141-006	3
8	Water Relay	123-000	1	22	Screw, 6-32 x 1/4"	001-133	4
9	Routing Valve	014-007	1	23	Needle Flow Valve	121-000	1
10	Plug, 1/2"	010-013	2	24	Automatic Handpiece Holder	122-049	1
11	Chassis	101-014	1	25	Screw, 4 x 3/4"	001-063	2
12	Plastic Elbow, 1/16"	022-039	1	26	Screw, 6-32 x 1/2"	001-072	4
13	Gauge	026-003	1	27	Washer, Flat #6	003-032	4
14	Nut, 6-32	002-004	2	28	Foot Control	126-000	1
15	Screw, 4 x 3/8"	001-135	4	29	Cart Frame	108-016	1
16	Faceplate Extension	101-024	2	30	Screw, 1/4-20 x 7/8"	001-011	4
17	Screw, 4 x 1/2"	001-063	4	31	Lock Washer, #6	003-033	4
18	Faceplate	101-015	1	32	Screw, 10-32 x 3/8"	001-070	1
19	Vacuum Holder	141-038	1	33	Handpiece Adjustment Block	113-021	1
20	3-Way Toggle Valve	120-008	2	34	Ferrule, 1/4"	117-008	1
21	Universal Handpiece Holder	141-006	3	35	Ferrule, 1/8"	117-007	1
22	Screw, 6-32 x 1/4"	001-133	4	36	Shuttle Valve	014-012	1
23	Needle Flow Valve	121-000	1	37	Plastic Tee, 1/8"	022-049	1
24	Automatic Handpiece Holder	122-049	1	38	Plastic Tee, 1/16"	022-033	1
25	Screw, 4 x 3/4"	001-063	2	39	Storage Organizer Assembly	110-016	1
26	Screw, 6-32 x 1/2"	001-072	4	40	Lockwasher, 1/2"	003-017	4
27	Washer, Flat #6	003-032	4	41	Syringe	113-400	1
28	Foot Control	126-000	1				

**AUTOMATIC HYGIENIST'S CART
A-2150**

Rev 2 4/4/10