



ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY
1	Chassis.....	101-014	1	10	Screw, 4 x 3/8".....	001-135	6	19	Plug, 10-32.....	001-001	1
2	Plug, 1/2".....	010-013	2	11	Universal Handpiece Holder.....	141-006	3	20	Shuttle Valve.....	014-012	1
3	Plastic Elbow, 1/16".....	022-039	1	12	3-way Toggle Valve.....	120-000	1	21	Plastic Tee, 1/8" x 1/16" x 1/8".....	022-049	1
4	Gauge.....	026-003	1	13	Ferrule, 1/4".....	117-008	1	22	Syringe Block.....	113-021	1
5	Nut, 6-32.....	002-004	2	14	Barb, 10-32 x 1/8".....	022-009	1	23	Water Relay.....	123-000	1
6	Washer, Flat, #6.....	003-032	4	15	Washer, nylon, 1/4".....	003-001	1	24	Screw, 10-32 x 3/8".....	001-070	1
7	Faceplate.....	101-015	1	16	4-way Cross.....	023-008	1	25	Needle Flow Valve.....	121-000	1
8	Screw, 6-32 x 1/4".....	001-133	2	17	Barb, 10-32 x 1/16".....	022-010	1	26	Unit Cover (optional).....	101-016	1
9	Routing Valve.....	014-006	1	18	Ferrule, 1/8".....	117-007	1	27	Screw, 6-32 x 3/8".....	001-134	2
								28	Syringe.....	113-400	1