



| ITEM | PART DESCRIPTION          | PART NUMBER | QTY | ITEM | PART DESCRIPTION                 | PART NUMBER | QTY | ITEM | PART DESCRIPTION                 | PART NUMBER | QTY |
|------|---------------------------|-------------|-----|------|----------------------------------|-------------|-----|------|----------------------------------|-------------|-----|
| 1    | Chassis .....             | 101-014     | 1   | 11   | Universal Handpiece Holder ..... | 141-006     | 2   | 21   | Plastic Tee, 1/16".....          | 022-033     | 1   |
| 2    | Plug, 1/2".....           | 010-013     | 2   | 12   | Automatic Handpiece Holder.....  | 122-049     | 1   | 22   | Syringe Block .....              | 113-021     | 1   |
| 3    | Plastic Elbow 1/16" ..... | 022-039     | 1   | 13   | Screw, 4" x 3/4".....            | 001-063     | 2   | 23   | Water Relay.....                 | 123-000     | 1   |
| 4    | Gauge .....               | 026-003     | 1   | 14   | Ferrule, 1/4".....               | 117-008     | 1   | 24   | Routing Valve.....               | 014-007     | 1   |
| 5    | Nut, 6-32 .....           | 002-004     | 2   | 15   | Barb, 10-32 x 1/8".....          | 022-010     | 1   | 25   | Needle Flow Valve.....           | 121-000     | 1   |
| 6    | Washer Flat. #6.....      | 003-032     | 4   | 16   | Plastic Washer, 1/4".....        | 003-001     | 1   | 26   | Screw, 10-32 x 3/8".....         | 001-070     | 1   |
| 7    | Faceplate .....           | 101-015     | 1   | 17   | Barb, 10-32 x 1/16".....         | 022-009     | 1   | 27   | Tray Mount.....                  | 105-074     | 1   |
| 8    | Screw, 6-32 x 1/4".....   | 001-133     | 4   | 18   | Ferrule, 1/8".....               | 117-007     | 1   | 28   | Tray, 9/3/4" x 13 1/2"/S/S ..... | 010-001     | 1   |
| 9    | 3-Way Toggle Valve.....   | 120-007     | 2   | 19   | Shuttle Valve.....               | 014-012     | 1   | 29   | Mount Pivot Control Arm .....    | 106-203     | 1   |
| 10   | Screw, 4" x 3/8".....     | 001-135     | 1   | 20   | Plastic Tee, 1/8" .....          | 022-049     | 1   | 30   | Syringe .....                    | 113-400     | 1   |

**AUTOMATIC HANDPIECE CONTROL  
USED WITH A-2500**