



ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY
1	Backplate.....	103-300	1	8	Plastic Elbow, 1/16".....	022-039	1	15	Plastic Washer, 1/4".....	003-001	
2	Handpiece Manifold.....	023-023	1	9	Gauge.....	026-002	1	16	Barb, 10-30 x 1/16".....	022-009	
3	Needle Flow Valve.....	121-000	1	10	Universal Handpiece Holder.....	141-006	1	17	Ferrule, 1/8".....	117-007	
4	3-way Toggle Valve.....	120-000	2	11	Automatic Handpiece Holder....	122-049	3	18	Automatic Handpiece Block.....	125-030	1
5	Faceplate.....	103-301	1	12	Syringe.....	113-400	6	19	Water Relay Valve.....	123-000	2
6	Screw, 6-32 x 1/4".....	001-133	4	13	Ferrule, 1/4".....	117-008					
7	Screw, #4 x 3/8".....	001-135	8	14	Barb, 10-30 x 1/8".....	022-010					

**AUTOMATIC HANDPIECE CONTROL
A-400**