



ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY
1	Cup Fill Spout.....	132-068	1	17	Bushing.....	013-017	1	33	Electric Tgl. Switch, Low Voltage.....	016-020	1
2	Body, Cuspidor Spouts, Nickel..	132-065	1	18	Screw, 10-32 x 3/8".....	001-070	4	34	Universal Handpiece Holder.....	141-023	1
3	Plastic Elbow, 1/8".....	022-038	4	19	Vacuum Cover.....	132-109	1	35	Autoclavable S/E Valve.....	132-102	1
4	Cup Filler Valve.....	132-018	2	20	Screw, 1/4" -20 x 1-1/4".....	001-169	1	36	Vacuum Handpiece Holder.....	141-037	1
5	Bowl Rinse Spout.....	132-067	1	21	Post Clamp.....	105-088	1	37	Autoclavable Standard HVE Valve...112-262	112-262	1
6	Bowl Rinse Valve.....	132-020	1	22	Post Mount Bracket.....	132-104	1	38	Drain Assembly.....	132-041	1
7	Screw, 6-32 x 5/8".....	001-159	2	23	Screw, 8-32 x 1/4".....	001-149	2	39	Holder Bar.....	132-111	1
8	Flat Washer 5/8 x 1.....	003-017	1	24	Mounting Bracket, CV Cannister....	131-039	1	40	Cuspidor Frame.....	132-102	2
9	Hex Nut, 5/8"-24.....	002-519	1	25	Strain Relief.....	132-108	2	41	Screw, 10-32 x 3/4".....	001-034	2
10	Air Reservoir.....	132-044	1	26	HVE Plug.....	010-009	1	42	Cover.....	132-101	1
11	Signal Valve.....	132-035	1	27	Vacuum Cannister (uses 131-052 screen)	131-021	1	43	O-Ring, Ceramic Bowl.....	017-150	1
12	Cup Holder Drain Barb.....	132-047	1	28	Washer, Flat 7/16".....	003-035	2	44	Bowl, Ceramic.....	132-110	1
13	Screw, 10-32 x 1/4".....	001-087	4	29	Washer, Flat, (.515"ID, .875"OD).....	003-059	1	45	Strainer.....	132-008	1
14	Plug, .652-1/2".....	010-046	1	30	1/4" Quick Disconnect.....	024-003	1				
15	Bleed Valve.....	132-122	1	31	Electric Wire, Low Voltage.....	016-008	1				
16	Drain Shaft.....	132-103	1	32	Washer, SAE Flat 7/16" S.....	003-047	1				

**ASEPSIS CUSPIDOR
C-3250**