



ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY	ITEM	PART DESCRIPTION	PART NUMBER	QTY
1	Hinge	107-083	2	12	Holder Bar Mount	107-007	1	23	Indicator	026-020	1
2	Hex Nut 6-32	002-004	8	13	Holder Bar	107-075	1	24	Momentary Valve	120-005	1
3	Cover	107-112	1	14	Automatic Holder	122-100	3	25	Shuttle Valve	014-012	1
4	Chassis	107-110	1	15	Syringe Holder	141-100	1	26	Plastic Tee, 1/16"	022-033	1
5	Automatic Handpiece Block	125-030	1	16	Screw, 8 x 1/2"	001-018	8	27	Screw, 8-32 x 3/8"	001-151	4
6	Handpiece Manifold	023-001	1	17	Screw, 10-32 x 3/8"	001-070	1	28	Syringe	113-400	1
7	Gauge	026-003	1	18	Ferrule, 1/4"	117-008	2	29	Screw, 10-32 x 1/2"	001-166	2
8	Needle Flow Valve	121-000	4	19	Barb, 10-32 x 1/8"	022-010	4	30	Screw, 1/4"-20 x 1"	001-103	1
9	3-way Toggle Valve	120-000	3	20	Plastic Washer, 1/4"	003-001	13	31	Pivot Adapter Bushing	105-163	1
10	Base	107-113	1	21	Barb, 10-32 x 1/16"	022-009	9	32	Pivot Spacer	105-164	2

**AUTOMATIC HANDPIECE CONTROL
SO-3354**

Rev 2 2/10/10