



AMERICAN MADE

BDS



CABINET MOUNTED SYSTEMS

- Flexible options
- Specialty solutions
- Cabinet mounted hand controls



Model
SD-3302

CABINET MOUNTED SYSTEMS

With an industry-leading assortment of positioning options, BDS specializes in cabinet mounted handpiece controls. Whether the requirement is rear delivery, side delivery, arm mounted, or shelf mounted, BDS has every option needed to provide specialty solutions for every practice and preference.

Model ST-4250

- Doctor's side = Three-handpiece automatic with air/water controls
- Assistant's side = Three-way syringe, autoclavable S/E and HVE valves
- System Purge feature flushes entire water system and replaces with clean water
- Mounts to underside of counter top with included hardware
- Easy left-right positioning



Model SC-4200

- Three-handpiece automatic with air/water controls
- System Purge feature flushes entire water system and replaces with clean water
- Mounts to underside of counter top with included hardware
- Side mount, recessed controls



Model A-5150

- Three-way syringe, autoclavable S/E and HVE valves
- Mounts to underside of counter top with included hardware



Model

PM-3102, RD-3152 *(shown)*

- Three-handpiece automatic with air/water controls
- Autoclavable S/E and HVE valves (SO-3354)
- System Purge feature flushes entire water system and replaces with clean water
- System extends 12" from recessed position and pivots 30° left-right
- Mounts to underside of counter top with included hardware



Model

SD-3302

- Three-handpiece automatic with air/water controls
- System Purge feature flushes entire water system and replaces with clean water
- Heavy-duty, pneumatic air-brake flex-arm (9" or 24" rigid and 24" flex section with 25" vertical travel)
- Side mount, recessed controls
- Mounts to side or rear of cabinet with included hardware



Model

A-230, A-250 *(shown)*

- Two-handpiece with air/water controls (one wet and one dry)
- Automatic handpiece selection (A-250)
- Designed for under-shelf mounting
- A-230 = Without automatic handpiece selection
- A-250 = Automatic handpiece selection



503-538-8756

info@beaverstatedental.com

FEATURES OF CABINET MOUNTED SYSTEMS	SC-4250	SC-4200	A-5150	A-230	A-250	SO-3304	SO-3354	SD-3302
2 handpiece automatic (1 wet, 1 dry)Non-retracting 2 2					■			
2 handpiece manual (1 wet, 1 dry)				■				
3 handpiece automatic	■	■				■	■	■
3-way air/water syringe	■	■	■	■	■	■	■	■
Master on/off toggle with indicator	■	■				■	■	■
Pressure adjustment for each handpiece	■	■		■	■	■	■	■
Easy to read pressure gauge	■	■		■	■	■	■	■
Non-retracting water coolant, adjustment & purge	■	■				■	■	■
Non-retracting water coolant & adjustment				■	■			
Air coolant adjustment	■	■				■	■	■
Asepsis handpiece tubing	■	■	■	■	■	■	■	■
Vacuum canister with removable screen (131-052)						■		
Vacuum canister with removable screen (131-053)	■		■					
1 Autoclavable HVE and saliva ejector valve w/ tubing	■	■	■				■	
Asepsis flex-arm system with pneumatic brake								■
2" dia. wall mount bracket								■
Slide-out mount with 15° pivot left and right						■	■	
Standard disc foot control	■	■		■	■	■	■	■
Master valve assembly with filter and regulators	■	■				■	■	■
7' umbilical tubing	■	■	■	■	■	■	■	■

Please contact BDS for technical cut-out templates.



www.bdsystems.us

115 S. Elliott Road, Newberg, OR 97132

800.237.2303

503.538.8756

info@beaverstatedental.com