



AMERICAN MADE

BDS



Cuspidor SYSTEMS

Time honored solution
Clean and cost effective





Model
C-4250



CUSPIDOR SYSTEMS

Tried and true, BDS provides a cuspidor option what will integrate into any operatory. It's clean, it's cost effective, and our cuspidor system will retrofit onto any 2" post. You want them, we have them.

Model
C-4200



Model
C--4250



503-538-8756

info@beaverstatedental.com

FEATURES OF CUSPIDOR SYSTEMS	C-4200	C-4250
Easy clean asepsis design	■	■
Ceramic bowl	■	■
Easily removed autoclavable cup filler and bowl rinse spout	■	■
Adjustable timed bowl rinse	■	■
Gravity drain with dual venting	■	■
Vacuum canister w/ removable screen (uses 131-052 screen)		■
1 Autoclavable HVE and S/E valve with tubing		■
1/4" Female water quick disconnect		■
Low voltage electric vacuum switch		■
2" dia. post mount	■	■
OPTIONS FOR CUSPIDORS		
Additional umbilical (per foot)	■	
Additional umbilical (per foot)		■
132-121 Timed cup filler	■	■
Silcryn tubing		■
3-way syringe with tubing and holder	■	■

Please contact BDS for technical cut-out templates.



www.bdsystems.us
 115 S. Elliott Road, Newberg, OR 97132
 800.237.2303

503-538-8756
 info@beaverstatedental.com