



AMERICAN MADE

BDS



HANDPIECE CONTROL SYSTEMS

Simple & compact controls
Cost effective solutions
Versatile for all practices & labs



Model
S-5325



Model
S-330

HANDPIECE CONTROL SYSTEMS

When an operator is down, but a completely new package is not needed, BDS has retrofitting options to get that operator back up and running quickly, efficiently, and economically. From single handpiece

Model S-330

- Three-handpiece automatic with air/water controls
- Rugged aluminum cover
- Side mount, recessed controls
- Non-retracting water coolant and purge



Unit head dimensions (w x d) 14" x 13"

Model S-5325

- Two-handpiece automatic (1 wet, 1 dry)
- Easy to read pressure gauge
- Vacuum canister with removable screen (131-053)
- Telescoping arm with breakaway feature, 29" to 40" horizontal reach
- Non-retracting water coolant



Unit head dimensions (w x d) (11" x 6½")

Model A-300

- Three-handpiece automatic with air/water controls
- Pressure adjustment for each handpiece
- Non-retracting water coolant and purge



operation to a three-handpiece control with an air/water syringe, BDS's wide range of simple and compact controls will fill the need of any dental practice or lab.

Model

A-220, A-222, A-200, A-202

- Two-handpiece automatic
- Pressure adjustment for each handpiece
- Non-retracting water coolant for 2 handpieces A-222
- Non-retracting water coolant for 1 handpiece (1 wet, 1 dry) A-220
- A-220 = 2 Handpiece automatic
- A-222 = 2 Handpiece automatic both wet (not shown)
- A-200 = 2 Handpiece manual
- A-202 = 2 Handpiece manual both wet (not shown)



Dimensions (w x h x d)
6 1/8" x 5 3/4" x 3 1/2"



Unit head dimensions (w x h x d)
7" x 5 3/4" x 3 1/2"

Model

A-400

- Three-handpiece automatic
- Pressure adjustment for each handpiece



Unit head dimensions (w x h x d)
6 1/8" x 5 3/4" x 3 1/2"

Model

A-100


- Single handpiece control with coolant
- Handpiece pressure adjustment
- Non-retracting water coolant



Model

A-101

- Single handpiece control with out coolant
- Handpiece pressure adjustment
- Easy to read pressure gauge



Model

A-050

- Single handpiece control with out coolant
- Handpiece pressure adjustment
- Easy to read pressure gauge

503-538-8756

info@beaverstatedental.com



www.bdsystems.us
115 S. Elliott Road, Newberg, OR 97132
800.237.2303

503-538-8756
info@beaverstatedental.com