




AMERICAN MADE

BDS



MOBILE SYSTEMS

Versatile and efficient
Wide range of available options
Robust reliability



Model
S-4350



Model
A-3150

MOBILE SYSTEMS

BDS is the go-to supplier for Mobile Carts. Mobile Systems offer inherent advantages, including positioning freedom, a wide range of available options, and robust reliability.

Model A-3150

- Doctor's side = Three-handpiece automatic with air/water syringe
 - Effortless left/right swing design
 - Easy to read controls
 - System purge feature flushes entire water system and replaces with clean water
 - Assistant's side = Three-way syringe, autoclavable S/E and HVE valves
 - 16" x 22" Star Board ST work surface
- A-3100 = No vacuum
 - A-3150 = With vacuum



Model A-3100

- Assistant's side only = Three-handpiece automatic with air/water syringe
 - Effortless left/right swing design
 - Easy to read controls
 - System purge feature flushes entire water system and replaces with clean water
 - 16" x 22" Star Board ST work surface
- A-3100 = No vacuum
 - A-3150 = With vacuum



Model

S-4100 and S-4150 *(shown)*

- Three-handpiece automatic with air/water syringe
- Autoclavable S/E and HVE valves (S-4150)
- Side mount, recessed controls
- System purge feature flushes entire water system and replaces with clean water

■ S-4100 = No vacuum (Inset)

■ S-4150 = With vacuum



Model

A-4500 and A-4550 *(shown)*

- Three-way syringe, autoclavable S/E and HVE valves (A-4550)
- Plastic storage organizer assembly
- 16" x 22" Star Board ST work surface

■ A-4500 = Cart only (Inset)

■ A-4550 = With vacuum and syringe



503-538-8756

info@beaverstatedental.com

MOBILE SYSTEMS

BDS is the go-to supplier for Mobile Carts. Mobile Systems offer inherent advantages, including positioning freedom, a wide range of available options, and robust reliability.



Model S-4250

- Doctor's side = Three-handpiece automatic with air/water syringe
- Side mount, recessed controls
- System purge feature flushes entire water system and replaces with clean water
- Assistant's side = Three-way syringe, autoclavable S/E and HVE valves
- Plastic storage organizer assembly
- 16" x 22" Star Board ST work surface

Model
A-3151

- Doctor's side = Three-handpiece automatic with air/water syringe
 - Easy to read controls
 - System purge feature flushes entire water system and replaces with clean water
 - Assistant's side = Three-way syringe, autoclavable S/E and HVE valves
 - 16" x 22" Star Board ST work surface
- A-3101 = No vacuum
 - A-3151 = With vacuum



Model
A-3101

- Doctor's side only = Three-handpiece automatic with air/water syringe
 - Easy to read controls
 - System purge feature flushes entire water system and replaces with clean water
 - 16" x 22" Star Board ST work surface
- A-3101 = No vacuum
 - A-3151 = With vacuum



503-538-8756

info@beaverstatedental.com

MOBILE SYSTEMS

BDS is the go-to supplier for Mobile Carts. Mobile Systems offer inherent advantages, including positioning freedom, a wide range of available options, and robust reliability.

Model A-3153

- Three-handpiece automatic with air/water syringe
- Autoclavable S/E and HVE valves
- 16" x 22" Star Board ST work surface
- Easy to read controls
- System purge feature flushes entire water system and replaces with clean water

■ A-3153 = With vacuum

■ A-3103 = No vacuum



Model A-3103

- Three-handpiece automatic with air/water syringe
- 16" x 22" Star Board ST work surface
- Easy to read controls
- System purge feature flushes entire water system and replaces with clean water

■ A-3103 = No vacuum

■ A-3153 = With vacuum



Model
S-4350

- Three-handpiece automatic with air/water syringe
- Autoclavable S/E and HVE valves
- 14" x 9" removable work surface
- Side mount, recessed controls
- System purge feature flushes entire water system and replaces with clean water

- S-4300 = No vacuum
- S-4350 = With vacuum



Model
S-4300

- Three-handpiece automatic with air/water syringe
- 14" x 9" removable work surface
- Side mount, recessed controls
- System purge feature flushes entire water system and replaces with clean water

- S-4300 = No vacuum
- S-4350 = With vacuum



503-538-8756

info@beaverstatedental.com

FEATURES OF MOBILE SYSTEMS	S-4100	S-4150	S-4250	S-4300	S-4350	A-3100	A-3101	A-3103	A-3150	A-3151	A-3153	A-4500	A-4550
3 handpiece automatic	■	■	■	■	■	■	■	■	■	■	■		
3-way air/water syringe	■	■	■	■	■	■	■	■	■	■	■		■
Master on/off toggle with indicator	■	■	■	■	■	■	■	■	■	■	■		
Pressure adjustment for each handpiece	■	■	■	■	■	■	■	■	■	■	■		
Easy to read pressure gauge	■	■	■	■	■	■	■	■	■	■	■		
Non-retracting water coolant adjustment and purge	■	■	■	■	■	■	■	■	■	■	■		
Asepsis handpiece tubing	■	■	■	■	■	■	■	■	■	■	■		
16" x 22" Star Board ST work surface dolphin gray			■			■	■	■	■	■	■	■	■
Removable Starboard™ tray Insert				■	■								
Painted H-frame with height adjustment, 27" to 35"	■	■	■	■	■	■	■	■	■	■	■	■	■
Storage organizer assembly			■									■	■
Vacuum canister with removable screen (131-052)									■	■	■		■
Vacuum canister with removable screen (131-053)		■	■		■								
Autoclavable HVE and saliva ejector valve with tubing		■	■		■				■	■	■		■
Standard disc foot control	■	■	■	■	■	■	■	■	■	■	■		
Swing for easy left/right conversion			■						■				

Please contact BDS for technical cut-out templates.



www.bdsystems.us
 115 S. Elliott Road, Newberg, OR 97132
 800.237.2303

503-538-8756
 info@beaverstatedental.com