




AMERICAN MADE

BDS



PANEL MOUNT SYSTEMS

Incredibly simple and efficient
Cost effective in a wide range of choices
Easy serviceability



Model
PM-4550

PANEL MOUNT SYSTEMS

Efficient, cost-effective, and vastly customizable, the advantages of panel mount delivery systems are endless. Easy to service, and built with BDS legendary simplicity and durability, the panel mount systems will keep the operatory, and entire office, running smoothly.

Model PM-4551

- Three-handpiece automatic with air/water controls
- Dual mounting options -
Panel Mount / Under Counter Mount
- Rotating dual holder bars for multi positioning
- Autoclavable SE and HVE



Panel mount dimensions (w x h x d)
9¾" x 2¾" x 9⅞"

Model PM-4550

- Three-handpiece automatic with air/water controls
- Rigid panel mount holder bar
- Autoclavable SE and HVE



Panel mount dimensions (w x h x d)
9¾" x 2¾" x 9⅞"

Model PM-4501

- Three-handpiece automatic with air/water controls
- Dual mounting options -
Panel mount / Under counter mount
- Rotating holder bar for multi positioning



Panel mount dimensions (w x h x d)
9¾" x 2¾" x 9⅞"

Model
PM-4500

- Three-handpiece automatic with air/water controls
- Rigid panel mount holder bar



Panel mount dimensions (w x h x d)
9¾" x 2¾" x 9⅞"

Model
PM-400

- Three-handpiece automatic
- Master ON/OFF toggle with indicator
- Water ON/OFF toggle with flow adjustment
- Easy read pressure gauge



Panel mount dimensions (w x h x d)
7½" x 2" x 5⅞"

Model
PM-200, PM-220 *(not shown)*

- Two-handpiece control (1 wet, 1 dry)
- Easy read pressure gauge
- Standard disc foot control
- PM-200 = Manual Control
- PM-220 = Automatic Control (not shown)



Panel mount dimensions (w x h x d)
6⅜" x 2" x 5⅞"

503-538-8756

info@beaverstatedental.com

FEATURES OF PANEL MOUNT SYSTEMS	PM-4551	PM-4550	PM-4501	PM-4500	PM-400	PM-220	PM-200
2 handpiece automatic						■	
3 handpiece automatic	■	■	■	■	■		
2 handpiece manual							■
3-way air/water syringe	■	■	■	■	■	■	■
Dual 3-way air/water syringes	■						
Master on/off toggle with indicator	■	■	■	■	■		
Pressure adjustment for each handpiece	■	■	■	■	■	■	■
Easy to read pressure gauge	■	■	■	■	■	■	■
Non-retracting water coolant adjustment, & purge	■	■	■	■			
Non-retracting water coolant with adjustment					■	■	■
Air coolant adjustment	■	■	■	■			
Vacuum canister with removable screen (131-052)	■	■					
Asepsis handpiece tubing	■	■	■	■	■	■	■
Standard disc foot control	■	■	■	■	■	■	■
7' Supply line							■
7' umbilical tubing	■	■	■	■	■	■	

Please contact BDS for technical cut-out templates.



www.bdsystems.us
 115 S. Elliott Road, Newberg, OR 97132
 800.237.2303

503-538-8756
 info@beaverstatedental.com