



AMERICAN MADE

BDS



NEW DESIGN

NEW POSSIBILITIES

** Compact and versatile*

** Narrow mounting footprint - under 12" in width*

** Cost effective solution*

** Multiple positioning preferences*



Model RD-3250

REAR CABINET MOUNTED SYSTEMS

If you are looking for the versatility of a rear cabinet mounted system, BDS offers a wide array of options for you. Configure it how you like and put it to work. BDS systems will keep your operatory working reliably, efficiently, and effectively.

Model

RD-3250 \$3,400.00

- Under counter mount duo system with left/right swing
- Three handpiece automatic with air/water controls
- Autoclavable SE and HVE
- Left/right capability
- Standard disc foot control
- Air coolant adjustment



Model

RD-3203 \$2,250.00

- Under counter mount doctor's system
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment



Model

RD-3201 \$2,400.00

- Under counter mount doctor's system
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment





Model

RD-3200 \$2,800.00

- Under counter mount doctor's system
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment



Model

RD-3253 \$2,600.00

- Under counter mount duo system
- Three handpiece automatic with air/water controls
- Autoclavable SE and HVE
- Standard disc foot control
- Air coolant adjustment



Model

RD-3251 \$2,700.00

- Under counter mount duo system
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment
- Autoclavable SE and HVE



503-538-8756

info@beaverstatedental.com

FEATURES OF REAR CABINET MOUNTED SYSTEMS	RD-3250	RD-3200	RD-3201	RD-3251	RD-3203	RD-3253				
3 handpiece automatic	■	■	■	■	■	■				
Non-retracting water coolant adjustment and purge	■	■	■	■	■	■				
Pressure adjustment for each handpiece	■	■	■	■	■	■				
Dual 3-way air/water syringes	■			■						
3-way air/water syringe		■	■		■	■				
Master on/off toggle with indicator	■	■	■	■	■	■				
Asepsis handpiece tubing	■	■	■	■	■	■				
Autoclavable, HVE saliva ejector valve & tubing	■			■		■				
Wall mounting plate										
Easy to read pressure gauge	■	■	■	■	■	■				
Standard disc foot control	■	■	■	■	■	■				
Vacuum canister with removable screen (131-052)	■			■		■				
Vacuum canister with removable screen (131-053)										
Master valve assembly with filter and regulators	■	■	■	■	■	■				
7' umbilical tubing	■	■	■	■	■	■				

Please contact BDS for technical templates.



www.bdsystems.us
 115 S. Elliott Road, Newberg, OR 97132
 800.237.2303

503-538-8756
 info@beaverstatedental.com