

# Cuspidor SYSTEMS

Time honored solution  
Clean and cost effective





Model  
C-4250



# CUSPIDOR SYSTEMS

Tried and true, BDS provides a cuspidor option what will integrate into any operator. It's clean, cost effective, and our cuspidor system will retrofit onto any 2" post. The simple, durable solution for patient evacuation.

Model  
**C-4200**



Model  
**C--4250**



800.237.2303

info@beaverstatedental.com



## Durability by Design

800.237.2303

info@beaverstatedental.com

FEATURES OF CUSPIDOR SYSTEMS	C-4200	C-4250
Easy clean asepsis design	■	■
Ceramic bowl	■	■
Easily removed autoclavable cup filler and bowl rinse spout	■	■
Adjustable timed bowl rinse	■	■
Gravity drain with dual venting	■	■
Vacuum canister w/ removable screen (uses 131-052 screen)		■
1 Autoclavable HVE and S/E valve with tubing		■
1/4" Female water quick disconnect		■
Low voltage electric vacuum switch		■
2" dia. post mount	■	■
OPTIONS FOR CUSPIDORS		
Additional umbilical (per foot)	■	
Additional umbilical (per foot)		■
132-121 Timed cup filler	■	■
Silcryn tubing		■
3-way syringe with tubing and holder	■	■

Please contact BDS for additional technical information.



[www.bdsystems.us](http://www.bdsystems.us)  
115 S. Elliott Road, Newberg, OR 97132  
800.237.2303

