

OVER THE PATIENT SYSTEMS

Clean and contemporary
Wide range of configurations
New operatory to replacement equipment



Model
S-3600

OVER-THE-PATIENT SYSTEMS

With a wide variety of arm-mounted systems, BDS has the Over-the-Patient solutions you need. From a full new operatory package with traditional three-handpiece

Model

S-3600, S-3606

- Three-handpiece automatic with air/water controls
- Recessed cover for procedure tray
- Ergonomic, easy-access controls
- Heavy-duty, pneumatic air-brake flex-arm (24" rigid and 24" flex sections with 24" vertical travel)
- System Purge feature flushes entire water system and replaces with clean water

■ S-3600 = Post mount (shown) ■ S-3606 = Top-post mount



Model

S-3630, S-3636

- Three-handpiece automatic with air/water controls
- Recessed cover for procedure tray
- Ergonomic, easy-access controls
- Heavy-duty, pneumatic air-brake flex-arm (24" rigid and 24" flex sections with 24" vertical travel)
- System Purge feature flushes entire water system and replaces with clean water
- Unit mount swivel tray 9-3/4" x 13-1/2"

■ S-3630 = Post mount (shown) ■ S-3636 = Top-post mount



Model

S-3420, S-3426

- Three-handpiece automatic with air/water controls
- Recessed cover for procedure tray
- Manual counter-balanced flex-arm (27.5" rigid and 26.5" flex sections with 27" vertical travel)
- Suitable for doctor or hygienist use

■ S-3420 = Post mount (shown) ■ S-3426 = Top-post mount



Unit Mount Touch Pad 107-145

- available on S-3600, S-3606, S-3630 & S-3636

system, to retro-fitting an existing chair with a European style system to simple hygiene controls, BDS gets operatories running quickly, efficiently, and reliably.

Model

S-6300, S-6307

- Three-handpiece automatic with air/water controls
- Ergonomic and functional Euro styling
- Recessed controls
- Heavy-duty, pneumatic air-brake flex-arm (24" rigid and 24" flex sections with 24" vertical travel)
- System Purge feature flushes entire water system and replaces with clean water

■ S-6300 = Post mount ■ S-6307 = Top-post mount (shown)



Model

S-3300, S-3306

- Three-handpiece automatic with air/water controls
- Rugged aluminum cover
- Side mount, recessed controls
- Heavy-duty, pneumatic air-brake flex-arm (24" rigid and 24" flex sections with 24" vertical travel)
- System Purge feature flushes entire water system and replaces with clean water

■ S-3300 = Post mount (shown) ■ S-3306 = Top-post mount



Model

S-2500, S-2506

- Two-handpiece automatic with air/water controls
- Manual counter-balanced flex-arm (27.5" rigid and 26.5" flex sections with 27" vertical travel)
- Recessed cover for procedure tray
- Designed for hygienist use

■ S-2500 = Post mount (shown) ■ S-2506 = Top-post mount



800.237.2303

info@beaverstatedental.com



Durability by Design

800.237.2303

info@beaverstatedental.com

FEATURES OF OVER-THE-PATIENT SYSTEMS	S-3600	S-3606	S-3630	S-3636	S-3300	S-3306	S-6300	S-6307	S-3420	S-3426	S-2500	S-2506
2 handpiece automatic control											■	■
3 handpiece automatic control	■	■	■	■	■	■	■	■	■	■		
3-way air/water syringe	■	■	■	■	■	■	■	■	■	■	■	■
Master on/off toggle	■	■	■	■	■	■	■	■	■	■	■	■
Pressure adjustment for each handpiece	■	■	■	■	■	■	■	■	■	■	■	■
Easy to read pressure gauge	■	■	■	■	■	■	■	■	■	■	■	■
Non-retracting water coolant with adjustment									■	■	■	■
Non-retracting water coolant adjustment and purge	■	■	■	■	■	■	■	■				
Air coolant adjustment	■	■	■	■	■	■	■	■				
Asepsis handpiece tubing	■	■	■	■	■	■	■	■	■	■	■	■
Asepsis flex-arm system with pneumatic brake	■	■	■	■	■	■	■	■				
Manual flex-arm system with friction brake									■	■	■	■
Unit mount post, 28"		■		■		■		■		■		■
Standard disc foot control	■	■	■	■	■	■	■	■	■	■	■	■
Utility center assembly with master valves & regulators	■	■	■	■	■	■	■	■	■	■	■	■
7' umbilical tubing	■	■	■	■	■	■	■	■	■	■	■	■

Please contact BDS for additional technical information.



www.bdsystems.us
115 S. Elliott Road, Newberg, OR 97132
800.237.2303

