



BDS



NEW DESIGN

NEW POSSIBILITIES

** Compact and versatile*

** Narrow mounting footprint - under 12" in width*

** Cost effective solution*

** Multiple positioning preferences*



Model RD-3250

| FEATURES OF REAR CABINET MOUNTED SYSTEMS | RD-3250 | RD-3200 | RD-3201 | RD-3251 | RD-3203 | RD-3253 |
|---------------------------------------------------|---------|---------|---------|---------|---------|---------|
| 3 handpiece automatic | ■ | ■ | ■ | ■ | ■ | ■ |
| Non-retracting water coolant adjustment and purge | ■ | ■ | ■ | ■ | ■ | ■ |
| Pressure adjustment for each handpiece | ■ | ■ | ■ | ■ | ■ | ■ |
| Dual 3-way air/water syringes | ■ | | | ■ | | |
| 3-way air/water syringe | | ■ | ■ | | ■ | ■ |
| Master on/off toggle with indicator | ■ | ■ | ■ | ■ | ■ | ■ |
| Asepsis handpiece tubing | ■ | ■ | ■ | ■ | ■ | ■ |
| Autoclavable, HVE saliva ejector valve & tubing | ■ | | | ■ | | ■ |
| Wall mounting plate | | | | | | |
| Easy to read pressure gauge | ■ | ■ | ■ | ■ | ■ | ■ |
| Standard disc foot control | ■ | ■ | ■ | ■ | ■ | ■ |
| Vacuum canister with removable screen (131-052) | ■ | | | ■ | | ■ |
| Vacuum canister with removable screen (131-053) | | | | | | |
| Master valve assembly with filter and regulators | ■ | ■ | ■ | ■ | ■ | ■ |
| 7' umbilical tubing | ■ | ■ | ■ | ■ | ■ | ■ |

Please contact BDS for technical templates.



www.bdsystems.us
 115 S. Elliott Road, Newberg, OR 97132
 800.237.2303

800.237.2303
 info@beaverstatedental.com