



Rear Cabinet Mounted Systems

At BDS our selection of Cabinet Mounted Systems offer the ideal solution. Compact, cost effective and versatile, our range of options can be configured to your personal positioning preferences.



- Cost Effective &
Versatile


Durability by Design
www.beaverstatedental.com
800.237.2303

Rear *Cabinet Mounted Systems*

Model

RD-4600 and RD-4650 (Shown)

- Under counter mount duo system (RD-4650)
- Under counter mount doctor's system (RD-4600)
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment
- Autoclavable SE and HVE (RD-4650)



Model

RD-4250 and RD-4251 (Shown)


- Wall/cabinet mount duo system (RD-4250)
- Wall/cabinet mount Assistant's system (RD-4251)
- Three handpiece automatic with air/water controls (RD-4250)
- Autoclavable SE and HVE
- Starboard work surface
- Industrial strength pneumatic support arm
- Left/right capability (RD-4250)



Model

RD-3103 and RD-3153 (Shown)

- Under counter mount duo system (RD-3153)
- Under counter mount doctor's system (RD-3103)
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment
- Autoclavable SE and HVE (RD-3153)



Model

RD-3102 and RD-3152 (Shown)

- Under counter mount duo system (RD-3152)
- Under counter mount doctor's system (RD-3102)
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment
- Autoclavable SE and HVE (RD-3152)



Model

RD-3100 and RD-3150 (Shown)

- Under counter mount duo system with left/right swing
- Three handpiece automatic with air/water controls
- Autoclavable SE and HVE
- Left/right capability
- Standard disc foot control
- Air coolant adjustment

RD-3100 = No vacuum

RD-3150 = With vacuum



Model

RD-3101 and RD-3151 (Shown)

- Under counter mount duo system (RD-3151)
- Under counter mount doctor's system (RD-3101)
- Three handpiece automatic with air/water controls
- Standard disc foot control
- Air coolant adjustment
- Autoclavable SE and HVE (RD-3151)



Features of Rear Cabinet Mounted Systems

	RD-3150	RD-3100	RD-3101	RD-3151	RD-3102	RD-3152	RD-3103	RD-3153	RD-4250	RD-4251	RD-4600	RD-4650
3 handpiece automatic	■	■	■	■	■	■	■	■	■		■	■
Non-retracting water coolant, adjustment and purge	■	■	■	■	■	■	■	■	■		■	■
Pressure adjustment for each handpiece	■	■	■	■	■	■	■	■	■		■	■
Dual 3-way air/water syringes	■			■		■			■			
3-way air/water syringe		■	■		■		■	■		■	■	■
Master on/off toggle with indicator	■	■	■	■	■	■	■	■	■		■	■
Asepsis handpiece tubing	■	■	■	■	■	■	■	■	■		■	■
Autoclavable, HVE saliva ejector valve & tubing	■			■		■		■	■	■		■
Wall mounting plate									■	■		
Easy to read pressure gauge	■	■	■	■	■	■	■	■	■		■	■
Standard disc foot control	■	■	■	■	■	■	■	■	■	■	■	■
Vacuum canister with removable screen (131-052)	■			■	■	■		■				■
Vacuum canister with removable screen (131-053)									■	■		
Master valve assembly with filter and regulators	■	■	■	■	■	■	■	■	■			
7' umbilical tubing	■	■	■	■	■	■	■	■	■	■	■	■



Hand Built
in the USA
7/2016

115 S. Elliott Rd.
Newberg, OR 97132
800.237.2303
www.beaverstatedental.com

