



Stools

BDS offers simple seating solutions for the dental operator. Both the Doctor's and Assistant's stool include height adjustment for proper orientation to the dental patient during a procedure. The Doctor's stool includes an adjustable backrest and the Assistant's stool has an adjustable torso support for proper positioning.

Model
DR-96



Model
AT-96




Durability by Design
www.beaverstatedental.com
800.237.2303

- Simply Seating

Features of Operatory Stools	DR-96	AT-96	DR-97	AT-97
Contoured seamless seat	■			■
Round seat, 15-1/2" x 4"		■		
Chrome plated back upright	■		■	
Contoured seamless semi-rectangular back upright	■			
Contoured seamless oval back			■	
Chrome plated, full padded body support		■		■
Chrome plated, pneumatic cylinder with height adjustment, 18"-22"	■		■	
Chrome plated, pneumatic cylinder with height adjustment, 20"-28"		■		■
5-leg, powder coated base	■	■	■	■
Twin wheel casters	■	■	■	■
Adjustable foot ring		■		■
Standard Nauga Soft upholstery	■	■	■	■
Contoured seamless seat with tilt			■	
Chrome plated back upright with tilt			■	
Ratcheting body support				■

Options of Operatory Stools	DR-96	AT-96	DR-97	AT-97
Ratcheting body support with plastic cover		■		
Hard floor casters, set of 5	■	■	■	■
Additional upholstery options may be available upon request	■	■	■	■
Ultraleather (in lieu of standard upholstery)	■	■	■	■
Taller cylinders available upon request	■	■	■	■



Hand Built
in the USA
3/2017

115 S. Elliott Rd.
Newberg, OR 97132
800.237.2303

www.beaverstatedental.com

